- ← ELEKTRODENPOSITIONERUNGSANLEITUNG
- ← ELECTRODE PLACEMENT GUIDE
- ← GUÍA PARA LA COLOCACIÓN DE ELECTRODOS
- ← GUIDE DE POSITIONNEMENT DES ÉLECTRODES
- ← GUIDA AL POSIZIONAMENTO DEGLI ELETTRODI
- ← OVERZICHT PLAATSING VAN DE ELEKTRODEN
- ← ELEKTRODEPLASSERINGSGUIDE
- ← ELEKTRODPLACERINGSGUIDE
- ← GUIA DE COLOCAÇÃO DE ELETRODO

	GEBRAUCHSANWEISUNG.....	3
	OPERATING INSTRUCTIONS	31
	MANUAL DE INSTRUCCIONES.....	59
	MODE D'EMPLOI	87
	ISTRUZIONI PER L'USO	115
	GEBRUIKERSHANDLEIDING.....	143
	BRUKSANVISNING	171
	BRUKSANVISNING	199
	INSTRUÇÕES DE FUNCIONAMENTO	227

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ADDITIONAL INFORMATION

8. ACCESSORIES	47
9. CARE INSTRUCTIONS.....	48
10. TROUBLESHOOTING.....	49
11. FREQUENTLY ASKED QUESTIONS.....	50
12. TECHNICAL DATA.....	51
KEY TO THE SYMBOLS.....	52
13. ELECTRODE PLACEMENT GUIDE.....	SEE COVER
14. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES.....	53

1. INTRODUCTION

CEFAR Medical AB has been a medical supplier of professional TENS Pain Relief products for more than 30 years. The company was established in Sweden in 1975 and is currently represented in large parts of the world.

The Cefar stimulators are widely used and endorsed by medical professionals and top athletes around the world.

TENS is effective, has no side effects and it's economical. Through clinical research, areas of application for TENS and EMS are rapidly expanding. Cefar is working actively to further develop the method towards a natural treatment alternative for the consumer.

More information about TENS (Transcutaneous Electrical Nerve Stimulation), EMS (Electrical Muscle Stimulation) and our products can be found on our web site: www.cefar.se

INTENDED USE:

CEFAR BASIC (1529) is a dual channel (four electrodes) nerve stimulator (TENS) with non-independent channels. CEFAR BASIC is designed for use by both professionals and home users. Its three preset pain programs provides relief of a wide range of pains, like back pain and neck pain. The stimulator has buttons for setting of program, amplitude and timer. An automatic key lock is activated 10 seconds after the amplitude has been set.

INTENDED USER:

The user of CEFAR BASIC can be a Health Care Professional or a patient. The device should be used indoors and may be used in a healthcare facility setting or in a home environment.

2. MEDICAL BACKGROUND

WHAT IS TENS?

TENS = Transcutaneous Electrical Nerve Stimulation. TENS is an effective method offering Pain Relief. TENS utilises the nervous system's own pain relief mechanisms to ease both acute and chronic pain.

WHEN IS TENS USED?

Good results are observed in acute and long-term pain conditions of many kinds, but above all when the pain originates in joints, the skeleton, muscles, skin, viscera, or nervous system. In some cases, TENS on its own may be sufficient, and in other cases it may provide a valuable complement to other forms of treatment.

HOW DOES IT WORK?

With the Pain Relief programs, impulses are sent through the nerves to block the pain impulses. Pain Relief is usually most efficient during stimulation, but the effect can also last after the program has finished. It also increases the circulation. The Pain Relief programs do not always treat the cause of pain. If the pain persists, seek medical advice.

WHEN AND HOW SHOULD I USE TENS?

The programs in the CEFAR BASIC can be used on any occasion when Pain Relief or Muscle Relaxation is needed. All the programs can be used as often as required. Guidelines for using the treatment programs with indication examples, suggested electrode placements and recommended treatment time are included in this manual.

3. PRECAUTIONARY MEASURES

3.1. CONTRAINDICATIONS

- Implanted electronic devices. Do not use the device on patients who have a cardiac pacemaker, implanted defibrillator, or other implanted electronic device, because this may cause electric shock, burns, electrical interference, or death. TENS for Undiagnosed Pain. Do not use the device as a TENS device on patients whose pain syndromes are undiagnosed.

3.2. WARNINGS

- Consult with physician. Consult with the patient's physician before using the device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- Skin condition. Apply stimulation only to normal, intact, clean, healthy skin.
- Long term effects. The long-term effects of chronic electrical stimulation are unknown.
- Stimulation location. Stimulation over Neck or Mouth. Do not apply stimulation over the patient's neck (especially the carotid sinus) or the patient's mouth, because this could cause severe muscle spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Stimulation across Chest. Do not apply stimulation across the patient's chest, because the introduction of electrical current into the chest may cause rhythm disturbances to the patient's heart, which could be lethal.
- Across the head. Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head, and electrodes should not be placed on opposite sides of the head.
- Stimulation over compromised skin. Do not apply stimulation over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins).
- Stimulation near cancerous lesions. Do not apply stimulation over, or in proximity to, cancerous lesions.
- Stimulation over eyes. Do not apply stimulation directly on the eyes.
- Environment. Electronic monitoring equipment. Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Bath or Shower. Do not apply stimulation when the patient is in the bath or shower. Do not apply stimulation in humid atmosphere exceeding 75% of relative humidity.
- Sleeping. Do not apply stimulation while the patient is sleeping.
- Driving or operating machinery. Do not apply stimulation while the patient is driving, operating machinery, or during any activity in which electrical stimulation or involuntary muscle contraction can put the patient at risk of injury.
- Electrosurgical equipment or defibrillators. Disconnect the stimulation electrodes before using electrosurgical equipment or defibrillators. Otherwise skin burns may be caused below the electrodes and the device might be destroyed.
- Magnetic Resonance Imaging. Do not wear electrode or the device during Magnetic Resonance Imaging (MRI) scans as this may result in metal overheating and causing skin burns in the area of the electrode.
- Flammable or explosive environment. Do not use the device in areas where there is a risk of fire or explosion, such as oxygen-rich environments, in the vicinity of flammable anaesthetics, etc.
- Power supply. Never connect stimulation cables to an external power supply as there is a risk of electric shock.
- Near other equipment. Do not use the device beside or stacked on top of any other equipment. If you must use it side by side or on top of another system, you should check that the device works properly in the chosen configuration.

- Miscellaneous. Electrodes for Single Patient. Do not share electrodes with other persons. All users
- should have individual set of electrodes to prevent undesirable skin reactions or disease transmission.
- Accessories. Use this device only with the leads, electrodes, and accessories recommended by the manufacturer. Use of other accessories may adversely affect the performance of the device or may result in stronger electromagnetic emissions or reduce the electromagnetic immunity of the device.
- No Modification. No modification of the equipment is allowed.

3.3. PRECAUTIONS

- Supervision. Use this device only under the continued supervision of a licensed practitioner. Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.
- Manufacturer. The manufacturer does not take any responsibility for any electrode placements other than recommended.
- Pregnancy. The safety of electrical stimulation during pregnancy has not been established.
- Skin irritation. Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (gel). The irritation may be reduced by using an alternate conductive medium or alternate electrode placement. Some patients may experience redness under the electrodes after a session. This redness usually disappears within a few hours. Advise the patient to consult the clinician if the skin redness does not disappear after a few hours. Do not start another stimulation session in the same area if the redness is still visible. Don't scratch the redness area.
- Heart Disease. Patients with suspected or diagnosed heart disease should follow precautions recommended by their physicians.
- Epilepsy. Patients with suspected or diagnosed epilepsy should follow precautions recommended by their physicians.
- Internal Bleeding. Use caution when the patient has a tendency to bleed internally, such as following an injury or fracture.
- After Surgery. Use caution following recent surgical procedures when stimulation may disrupt the patient's healing process.
- Over uterus. Do not place electrodes directly over the uterus or connect pairs of electrodes across the abdomen if you are pregnant. The reason is that, theoretically, the current could affect the foetus's heart (although there are no reports of it being harmful).
- Lack of sensation. Use caution if stimulation is applied over areas of skin that lack normal sensation. Don't apply stimulation on patient unable to express themselves.
- Stimulation. The stimulator should only be used with skin electrodes intended for nerve and muscle stimulation. Muscle soreness may occur after stimulation but usually disappears within a week.
- Hot casing or batteries. Under extreme use conditions, some parts of the casing might reach up to 109 °F (43 °C). Use caution when manipulating the batteries right after device use or when holding the device. There is no particular health risk associated with this temperature besides your comfort.
- Children. Keep this device out of the reach of children.
- Electrode Size. Do not use electrodes with an active area less than 16 cm², as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2mA/ cm².
- Strangulation. Do not wrap leadwires around your neck, and keep them out of the reach of children. Strangulation may result from entanglement in the leadwires.
- Tripping. Care should be used to avoid tripping on lead wires.
- Damaged Device or Accessories. Never use the device or any of its accessories if it is damaged (case, cables, etc.) or if the battery compartment is open as there is a risk of electric shock. Carefully inspect the lead wires and connectors prior to each use.

- Inspect Electrodes. Inspect electrodes before each use. Replace electrodes when they begin to deteriorate or lose adhesion. Poor contact between the electrodes and the patient's skin increases the risk of skin irritation or burns. Electrodes will last longer if used and stored according to instructions on electrode packaging. Attach the electrodes in such a way that their entire surface is in contact with the skin.
- Foreign bodies. Do not allow any foreign bodies (soil, water, metal, etc.) to penetrate the device and the battery compartment.
- Batteries. Do not carry batteries in a pocket, purse, or any other place where the terminals could become short-circuited (e.g. by way of paper clip). Intense heat could be generated and injury may result. Never open the battery cover during stimulation in order to avoid electrical shock. Remove the batteries from the device if you do not intend to use it for a prolonged period of time (more than 3 months). Extended storage of the batteries in the device might lead to batteries and device damage.
- Cable. The cable is best preserved if left attached to the stimulator between sessions. Do not jerk the cable or connection.
- Heat and Cold products. The use of heat or cold producing devices (e.g. electric heating blankets, heating pads or ice packs) may impair performance of the electrode or alter the patient's circulation/sensitivity and increase the risk of injury to the patient.
- Pulled muscles. Do not apply electrodes over pulled muscles. Using the stimulator on a previously extended muscle might further pull such muscle. The higher the stimulation intensity, the higher the risk to further overextend such muscle.
- Additional Precautions for TENS
 - TENS is not effective for pain of central origin, including headache.
 - TENS is not a substitute for pain medications and other pain management therapies.
 - TENS devices have no curative value.
 - TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- Effectiveness of TENS is highly dependent upon patient selection by a practitioner qualified in the management of pain patients.

3.4. DANGERS

- Electrodes. Any Electrode with a minimum active area of 16 cm² may be used with this device. Use of an electrode with an area less than 16 cm² can cause burns when the unit is used at higher intensities. Consult your clinician prior to using any electrode less than 16 cm².


3.5. ADVERSE REACTIONS

- Patients may experience skin irritation and burns beneath the stimulation electrodes applied to the skin.
- Patients may experience headache and other painful sensations during or following the application of electrical stimulation near the eyes and to the head and face.
- Patients should stop using the device and should consult with their physicians if they experience adverse reactions from the device.
- Precaution: Do not disconnect any stimulation cables during a session while the stimulator is switched on. Switch the stimulator off first. Always turn off the stimulator before moving or removing any electrodes during a session.
- Precaution: Do not use electrodes with an active area less than 16 cm², as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2mA/cm².
- Precaution: Do not apply stimulation in the vicinity of metal. Remove jewelry, body piercings, buckles or any other removable metallic product or device in the area of stimulation. Never use the electrodes contra-laterally, i.e. do not use two pins connected to the same channel on opposite segments of the body.

- Precaution: Never carry out an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, people of a nervous disposition may experience a vasovagal reaction. This reaction is connected with fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. A vasovagal reaction causes heart to slow down and blood pressure to drop, which can make you feel weak and faint. If this does occur, stop the stimulation and lie down with the legs raised until the feeling of weakness disappears (5 to 10 minutes).
- Precaution: Sudden temperature changes can cause condensation to build up inside the stimulator. To prevent this, allow it to reach room temperature before use.

4. OVERVIEW/CONTROL BUTTONS

1. ON/OFF

- Turns the stimulator on and off.
- Turns the stimulator off even when the key lock is activated.
- Can be used for terminating the stimulation at all times.

2. INCREASE (left and right channel)

- Increases the amplitude (intensity of stimulation). Press and hold the button to increase the amplitude continuously.

Note! Always increase the amplitude cautiously.

- Terminates AUTO stimulation mode.
- Increases the number of minutes when setting the timer (right button).
- Used when turning the program lock on and off (left button). (PROFESSIONAL USE)

3. DECREASE (left and right channel)

- Decreases the amplitude (intensity of stimulation). Press and hold the button to decrease the amplitude continuously.
- Deactivates the key lock.
- Terminates AUTO stimulation mode.
- Decreases the number of minutes when setting the timer (right button).
- Used when turning the program lock on and off (left button). (PROFESSIONAL USE)

4. AUTO (left and right channel)

- Starts the AUTO test mode.
- Confirms the selected amplitude level in AUTO test mode.

5. PROGRAM

- Selects programs 1-3.
- Pauses an ongoing program.
- Used when activating/deactivating the program lock.

6. TIMER

- Initiates the timer setting.


5. OPERATION


STEP-BY-STEP USE

1. INSERT THE BATTERIES

Insert the batteries (see chapter REPLACEMENT OF BATTERIES).


2. ATTACH THE ELECTRODES

A. Connect the electrodes to the cable.

B. Attach the electrodes to your body.


C. Connect the cable to the CEFAR BASIC.

3. SWITCH THE STIMULATOR ON

Press the **ON/OFF** button (⏻). This button can be used for terminating the stimulation at all times, even when the key lock is activated.

Always switch the stimulation off before removing electrodes from the skin.


4. SELECT A PROGRAM (1-3)

Press the **PROGRAM** button (P) repeatedly until the program of your choice is shown on the display.

Note! When selecting programs the amplitude must be set to 0.0 mA for both channels.

For further information on the programs, see section **PROGRAM GUIDE**.


6. REPLACEMENT OF BATTERIES

A battery symbol is shown on the display when the batteries are almost empty . As long as the stimulator is working normally you can continue the treatment. When stimulation feels weaker than usual or the stimulator turns off, it is time to replace the batteries.

If the stimulator is not used for some time (approximately 3 months), the batteries should be removed from the stimulator.

The stimulator operates on either two non-rechargeable 1.5 V AA batteries or two rechargeable 1.2 V AA batteries, recharged in a separate charger.

The endurance of the batteries depends on the amplitude and number of channels activated during a program being used. For example, if you use one channel at P1 with an amplitude of 30 mA the non-rechargeable batteries last approximately 11 hours and rechargeable 7 hours.

REPLACING THE BATTERIES

1. Turn off the stimulator.
2. Locate the battery compartment on the back of the stimulator.
3. Remove the battery cover by sliding it downwards about 1 cm and then lifting it upwards.
4. Remove the batteries.
5. Insert the new batteries correctly according to the polarity markings (+ and -) inside the battery compartment, see picture.
6. Replace the battery cover.
7. Dispose of the exhausted batteries in accordance with local and national regulations.

NOTE! Non-rechargeable batteries may not be charged due to the risk of explosion.

7. PROGRAM GUIDE

Each program is described with the indications it is best suited to treat, how the electrodes should be placed, and how long the treatment should last. To each example of indication is at least one Electrode Placement number, which refers to the suggested electrode placements shown in the cover of the manual.

No.	Program name	Frequency	Pulse duration	Treatment time	Preset time
1	Conventional TENS	80 Hz	180 μ s	No limit (minimum 30 minutes).	30 minutes
2	Burst TENS	2 Hz	180 μ s	20–45 minutes, 3 times a day. Remember that Burst TENS can cause muscle soreness.	30 minutes
3	Modulated pulse duration stimulation	80 Hz	70–180 μ s	No limit (minimum 30 minutes).	30 minutes

Program description	Examples of indications/ Electrodeplacement no.
<p>Conventional TENS (high-frequency stimulation) is the first choice for both acute and long-term pain, both neurogenic and nociceptive. Conventional TENS is based on the Gate Control theory, which states that electric stimulation of A-beta fibers inhibits impulse transfer in the pain pathways.</p> <p>As a rule, the electrodes should be placed on or near the painful area, or over an area segmentally related to the painful area. Adjust the amplitude so that the stimulation gives strong, but pleasant paraesthesia – tingling. When using high-frequency stimulation, it's important to ensure that the patient has normal sensory of touch in the area where the electrodes are placed.</p>	<p>Neck pain.....1 Shoulder pain2 Elbow pain6 Rheumatic pain.....7 Fractured rib12 Lumbago.....10,11 Menstrual pain.....13 Phantom limb pain.....14 Hip pain.....16 Osteoarthritic pain in the knee....17,18 Wound healing.....19</p>
<p>Burst TENS (low-frequency stimulation) is usually most effective for radiating (projected) pain in the arms and legs (rhizopathy), for conditions with reduced or changed sensory of touch, for deep muscular pain, or when the post-treatment effect of Conventional TENS is too short. Burst TENS treatment alleviates pain by stimulating muscles to release the body's own morphine-like substances, endorphins.</p> <p>Place the electrodes on a muscle in the painful area so that a visible contraction occurs, or on acupuncture points in the painful area. The stimulation should feel pleasant and give visible muscle contractions. Remember that the patient often feels the stimulation clearly before contractions become visible.</p>	<p>Mononeuropathy8 Central pain.....4, 5 Cervical rhizopathy.....3 Sciatica15 Knee pain.....18</p>
<p>Modulated pulse duration stimulation is a type of high-frequency stimulation where the pulse duration varies continuously. This can cause an undulating sensation, which may be more pleasant than a constant pulse duration. Use program 3 for pain alleviation and a massage effect on muscles like the trapezius.</p>	<p>Trapezius pain.....9 Lumbago.....10, 11 For more examples, see Program 1.</p>


8. ACCESSORIES

The electrodes will eventually wear out and need to be replaced. It is recommended to replace the electrodes after approximately 20-40 times of usage.

For purchase information, contact your Cefar dealer or visit www.cefarse.com

CEFAR BASIC comes with a neck strap and a belt clip, allowing you to wear the stimulator around your neck or on your belt and have your hands free during treatment.

The cables are best preserved if left attached to the stimulator between sessions.

9. CARE INSTRUCTIONS

Taking care of and cleaning the Cefar equipment is simple with the following instructions:

- Keep stimulator and accessories in the original case when they are not in use. It may, however, be practical to allow the electrodes to remain on the body between treatments. Carbon rubber electrodes can generally remain for 2–3 hours without the electrode gel drying out (does not apply to adhesive gel). They must then be taken off, washed, and dried before being applied again. This is especially important for persons with sensitive skin. In connection with stimulation, make sure that the electrodes are firmly in place.
- When using carbon rubber electrodes, use plenty of electrode gel and avoid drying out by applying tape around all the edges of the electrodes. Rinse the carbon rubber electrodes and the skin with water after use. Do not use detergent for the electrodes.
- Self-adhesive multi-use electrodes are re-moistened if necessary with a few drops of water and kept air-tight (in a plastic bag) on protective paper when they are not in use.
- Never expose the stimulator to water. Wipe it off with a damp cloth if necessary.
- Do not jerk cables or connections.
- The cables are best preserved if left attached to the stimulator between sessions.
- No other maintenance of the device is required. Service life of the device can vary depending on usage conditions. Typical service life is 7 years.
- Never service the device while in use.
- The device should be operated in temperatures between 10°C and 40°C, atmospheric pressures between 50 and 106 kPa, and relative humidity between 30% and 75%.
- The device should be transported and stored in temperatures between -40 °C and 70 °C, atmospheric pressures between 50 and 106 kPa and relative humidity between 10% and 90%.
- Contact manufacturer for assistance in setting up, using or maintaining the equipment or report events.

10. TROUBLESHOOTING

THE STIMULATION DOES NOT FEEL THE SAME AS USUAL

- Check that all settings are correct (see section STEP-BY-STEP-USE) and make sure that the electrodes are correctly placed.
- Slightly change the position of the electrodes.

THE STIMULATION FEELS UNPLEASANT

- The skin is irritated. For advice on skin care, see chapter PRECAUTIONARY MEASURE.
- The electrodes begin to lose their stickiness and do not stick properly to the skin. Moisten the adhesive surface with a few drops of water before placing on the skin.
- The electrodes are worn out and need to be replaced.
- There is insufficient electrode gel on the carbon rubber electrodes.
- Slightly change the position of the electrodes.

THE STIMULATION FEELS WEAK OR NOT AT ALL

- Check if the batteries need to be replaced, see chapter REPLACEMENT OF BATTERIES.
- Electrodes are too old and need replacement.

THE BROKEN CIRCUIT SYMBOL IS SHOWN ON THE DISPLAY SYMBOL

The broken circuit symbol indicates that the resistance is too high, or that a cable is broken.

- A too high resistance can be caused by a bad connection between the electrodes and your skin, or that the electrodes need to be replaced.
- A cable breakage can be checked by pressing the cable's pins against one another while increasing the amplitude for the corresponding channel to 11 mA. If the amplitude now drops to 0.0 mA and ⚡ starts flashing, the cable needs to be replaced.

Note! Never increase the amplitude above 20 mA when you check for cable breaks, since this can damage the stimulator.

THE STIMULATOR IS NOT WORKING

If the error symbol appears on the display when you start the stimulator, it means that the stimulator is broken and needs to be replaced.

Note! Do not use the stimulator – contact your Cefar dealer.

Cefar will only be responsible for service and repairs performed by Cefar or a distributor appointed by Cefar.

11. FREQUENTLY ASKED QUESTIONS (FAQ)

CAN ANYONE USE ELECTRICAL STIMULATION?

People with implanted electrical equipment for example a pacemaker and an intracardiac defibrillator must not be treated with electrical stimulation. Pregnant women should not use electrical stimulation during the first 12 weeks of the pregnancy. Read the safety precautions in this manual (PRECAUTIONARY MEASURES).

FOR HOW LONG TIME CAN I STIMULATE?

Follow the recommendations given by your health care provider. Generally high-frequency TENS (80 Hz) can be used without an upper time limit, but at least 30 minutes at each occasion. Low-frequency TENS (2 Hz) can cause sore muscles but normally 20-45 minutes 3 times a day is recommended.

CAN I USE TENS ON ALL KINDS OF PAIN?

TENS can be used for both acute and long time pain. General instructions can be read in chapter MEDICAL BACKGROUND. Contact your health care provider if the cause for and area of treatment is changed.

WHICH DISTANCE SHOULD I HAVE BETWEEN THE ELECTRODES?

It is recommended to have a distance of 3 to 30 cm between the electrodes.

FOR HOW LONG LAST THE ELECTRODES?

The self-adhesive electrodes last for approximately 20 to 40 occasions. The durability depends on how good the care and maintenance instructions are followed. Read the care and maintenance instructions delivered with the electrodes.

ARE THERE ANY SIDE EFFECTS WITH ELECTRICAL STIMULATION?

No, on the contrary to for example analgesic medicine, TENS stimulation is completely without side effects.

ARE THERE ANY PROOFS OF TENS WORKING AS AN ALLEVIATION OF PAIN?

TENS is an established and well-documented method that has been used by the medical service for many years. Research and medical studies are carried out all the time. Read more about current medical studies on Cefar's home page www.cefar.se.

12. TECHNICAL DATA

CEFAR BASIC is a nerve stimulator (TENS) with two non-independent channels and three preset programs.

Treatment with electrical stimulation requires the stimulation current to penetrate the resistance of the skin and the electrode, about 1500 ohms. CEFAR BASIC can penetrate this resistance and maintain a current of up to 60 mA. With a change in load from 100 to 1500 ohms, the stimulation current changes less than 10% from the set value.

The stimulator operates on either two non-rechargeable 1.5 V AA batteries or two rechargeable 1.2 V AA batteries, recharged in a separate charger.

CEFAR BASIC

Number of channels	2 (non-independent)
Constant current	Up to a resistance of 1500 ohm (increased load can reduce the maximum current)
Stimulation current/channel	0–60 mA (charge maximum: 10.8µC)
Waveform	Asymmetrical biphasic pulse, 100% compensated
Number of preset programs	3
Stimulation forms	Conventional Burst Modulated pulse duration
Max pulse duration	180 µs
Max frequency	80 Hz
Timer	0 to 60 min
Environment for storage, use and shipping	Temperature 10°C–40°C Air humidity 30%–75% Air pressure 700 hPa–1060 hPa
Power source	2 x 1.5 V AA non-rechargeable or 2 x 1.2 V AA rechargeable batteries
Current consumption for one channel, 80 Hz, 30 mA	150 mA
I r.m.s. max/channel	11 mA
Size	120 x 50 x 30 mm
Weight	ca. 180 g

14. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The CEFAR BASIC Electrotherapy System is intended for use in the electromagnetic environment defined below. The customer or the user of the CEFAR BASIC Electrotherapy System should assure that it is used in such an environment		
Emissions Tests	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The CEFAR BASIC Electrotherapy System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	CEFAR BASIC Electrotherapy System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Not Applicable - Battery powered	
Voltage fluctuations IEC 61000-3-3	Not Applicable - Battery powered	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CEFAR BASIC Electrotherapy System is intended for use in the electromagnetic environment specified below. The customer or the user of the CEFAR BASIC Electrotherapy System should assure that it is used in such an environment

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Risk assessment on the CEFAR BASIC Electrotherapy System indicates the compliance levels claimed are acceptable when ESD-precautionary measures are taken.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	Not Applicable - Battery powered Not Applicable - signal lines less than 3 meters	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+ 1kV differential mode (line to line) + 2kV common mode (line to ground)	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CEFAR BASIC Electrotherapy System requires continued operation during power mains interruptions, it is recommended that the CEFAR BASIC Electrotherapy System be powered from an uninterrupted power supply.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the AC mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CEFAR BASIC Electrotherapy System is intended for use in the electromagnetic environment defined below. The customer or the user of the CEFAR BASIC Electrotherapy System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	[V _i] V, where V ₁ = 3V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the CEFAR BASIC Electrotherapy System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \frac{3.5 \sqrt{P}}{V_1}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	[E ₁] V/m, where E ₁ = 3V/m	<p>$d = \frac{3.5 \sqrt{P}}{E_1}$ 80 MHz to 800 MHz</p> <p>$d = \frac{7 \sqrt{P}}{E_1}$ 800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CEFAR BASIC Electrotherapy System is used exceeds the applicable RF compliance level above, the CEFAR BASIC Electrotherapy System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CEFAR BASIC Electrotherapy System.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V_i] V/m.

Recommended separation distances between portable and mobile RF communications equipment and the CEFAR BASIC Electrotherapy System			
The CEFAR BASIC Electrotherapy System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CEFAR BASIC Electrotherapy System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CEFAR BASIC Electrotherapy System as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter P (W)	Separation distance according to frequency of transmitter d (m)		
	150 kHz to 80 MHz $d = \frac{3,5}{[V]} \sqrt{P}$ (where V ₁ = 3V)	80 MHz to 800 MHz $d = \frac{3,5}{[E_1]} \sqrt{P}$ (where E ₁ = 3V/m)	800 MHz to 2,5 GHz $d = \frac{7}{[E_1]} \sqrt{P}$ (where E ₁ = 3V/m)
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

FCC REQUIREMENTS

Part 15 of the FCC Requirements	
This device complies with Part 15 of the FCC Rules. Operation is subject to the following 2 conditions:	<ul style="list-style-type: none"> • This device may not cause harmful interference • This device must accept any interference received, including the interference that may cause undesired operation.
FCC ID	T9J-RN42
Contains Transmitter Module IC	6514A-RN42

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Photo © Studio Ekwall, Helsingborg AB